




AiDashboard™

User Guide

Ai-CURE-TECHNOLOGIES



**Ai-CURE TECHNOLOGIES**

URL: [www.myclinicaltrialdashboard.com](http://www.myclinicaltrialdashboard.com)


Call Technical Support: +1 (800) 675-4960 [7am to 7pm PST]

Email Technical Support: support@aicuretechnologies.com

Name:

Username:

Password:




Ai-CURE TECHNOLOGIES

- 1 About AiDashboard
- 2 Using the Events Section
- 3 Using the Patients Section
- 4 Using the Patient Chart Section
- 5 Analytics and Settings



## About the AiDashboard™

The AiDashboard™ allows you to track study participant performance. The dashboard provides real time adherence data for study participants, indicates which participants may need help, and suggests where you can best spend your time. Complete audit trails for each participant are also available.

# Quick Reference





# Training Study Participants



**IMPORTANT NOTE:** Please ask the study participant to plug their device into an electric power source and lock their device securely as soon as they return home.

# Participant & Study Coordinator Daily Use





Ai-CURE TECHNOLOGIES


- 1 About AiDashboard
- 2 Using the Events Section
- 3 Using the Patients Section
- 4 Using the Patient Chart Section
- 5 Analytics and Settings

# AiDashboard Login

Each study coordinator will be provided with a unique username and password to access the secure AiDashboard. To logon, please go to:

[www.myclinicaltrialdashboard.com](http://www.myclinicaltrialdashboard.com)

If you forget your password, please click on the orange “Forgot your password?” link, enter your email address and username on file, and a new password will be sent to your address.



# Events Section

When you first login, you will arrive at the EVENTS page which reports patient behavior such as missed, late, or skipped doses. Trial information (site, study drug, and participant group) and PI are displayed in the top header of the display. This EVENTS HOME PAGE is divided into three main panels:

## ACTION ITEMS

display items that require your attention, sorted by date.

The screenshot shows a list of patient events under the 'Events' tab. The events are listed by date, from today to the previous week. Each event is preceded by a small red icon indicating the nature of the event (e.g., missed dose, late dose). The list includes:

- Patient 8 - 8 has missed the previous 3 doses (Oct 4 2012)
- Patient 1 - 1 missed the 1st dose (Oct 3 2012)
- Patient 8 - 8 missed the 1st dose (Oct 2 2012)
- Patient 8 - 8 missed the 1st dose (Oct 1 2012)
- Patient 8 - 8 missed the 1st dose (Sep 30 2012)
- Patient 5 - 5 missed 1 pill of the 1st dose (Sep 30 2012)
- Patient 7 - 7 was 2 hours late for the 1st dose and took 1 pill of the dose manually (Sep 28 2012)
- Patient 3 - 3 manually took 4 mill of the 1st dose (Sep 28 2012)

PENDING

search patients by name

## PAPER STACK

displays more detailed information for a selected action item.

The screenshot shows a detailed view of a patient's event. It displays the date (30 Sunday September 2012) and the patient's name (Patient 5 - 5). The event is described as "Patient 5 - 5 missed 1 pill of the 1st dose." Below this, there is a section titled "Steps to be completed:" with several checkboxes. One checkbox is checked, and another is unchecked. There are also three small images showing different stages of taking a pill: "Pill identification" (pill in hand), "Pill in mouth" (pill in mouth), and "Swallow and complete" (pill swallowed). A "Watch Video" button is also present. At the bottom, there are two more checkboxes: one checked and one unchecked.

## COMPLETED TASKS

log tasks that have been completed in the dashboard.

The screenshot shows a list of completed tasks. Each task is listed with a date and a checkmark indicating it was completed. The tasks include:

- Completed steps for 8 - 8 from 2012-10-03
- Reviewed pill 3 for 3 - 3 from 2012-09-28
- Contacted 8 - 8 about 2012-10-03
- Added note for 8 - 8 about 2012-10-03
- Reviewed face identification for 7 - 7 from 2012-09-30
- Reviewed pill 3 for 5 - 5 from 2012-09-30

Thursday Oct 4 2012

Tuesday Oct 2 2012

Monday Oct 1 2012

Events Tasks Patients Date

Completed tasks

Settings

# Events Section - Action Items

The action items are reported in chronological order under the DATE tab, and provide information about events that require your attention, such as if a participant is late or has skipped or missed a dose. Once you open and begin addressing some of the tasks related to the ACTION ITEM, it will be moved to the PENDING list in the ACTION ITEMS panel.

The screenshot shows the Allure Technologies software interface. At the top, there is a navigation bar with the text "Pharmaceutical 2012 Trial | Site NY-1 | Med A | Group 1" and "PI: Dr. John Smith, PhD". To the right of the navigation bar are icons for "Events", "Tasks", "Patients", and "Date". The "Events" tab is currently selected. In the center of the screen, there is a calendar for September 2012, specifically focusing on Sunday, September 30, 2012. This date is highlighted in red. Below the calendar, there is a message: "Patient 5 - 5 missed 1 pill of the 1st dose." A detailed description follows: "Patient 5 - 5 skipped the 1 pill for their 1st dose. The patient started AlView and then missed taking some of their medication. You may want to check in with the patient to see if there is a problem." Below this message, there is a section titled "Steps to be completed:" which lists several tasks with checkboxes. One task is checked: "Patient 5 - 5 was correctly identified by the system." Another task is partially completed: "Pill 1 of 3 was taken correctly". Under this task, there is a sub-section titled "Pill 2 of 3 was taken correctly" with a checkbox that is also checked. There are four small images below this sub-section: "Pill identification" (showing a green outline), "Pill in mouth" (showing a green outline), "Swallow and complete" (showing a green outline), and "Watch Video" (showing a play button icon). Further down, there is another task: "Pill 3 of 3 produced errors. Please review." and a task at the bottom: "Contact Patient 5 - 5 to discuss Incident." On the left side of the screen, there is a sidebar titled "Events" which lists various patient actions in chronological order. Some entries are highlighted in pink, while others are in grey. The sidebar includes dates from October 1, 2012, to October 4, 2012. At the bottom of the sidebar, there is a search bar with the placeholder text "search patients by name" and a "Settings" link. On the right side of the screen, there is a panel titled "Completed tasks" which lists several tasks with checkmarks and dates from October 1, 2012, to October 4, 2012.

# Events Section - Paper Stack

When an ACTION ITEM is selected, the PAPER STACK will scroll down to reveal detailed information about an ACTION ITEM. The user may be asked to complete a series of tasks for each ACTION ITEM and check the box upon completion. When all tasks have been completed, the user may CLOSE EVENT, and the ACTION ITEM will collapse in the PAPER STACK and be logged. For each pill taken, de-identified video and images may be available. Clicking the icon for each pill will expand to show these images and video. You may click CONTACT PATIENT and send messages directly to a patient's device and / or write A NOTE (left of each step), to record information or patient communications related to a specific event.



The screenshot shows the Amurex Technologies software interface. At the top, there is a navigation bar with the text "Pharmaceutical 2012 Trial | Site NY-1 | Med A | Group 1" and "Pt: Dr. John Smith, PhD". To the right of the navigation bar are icons for "Events", "Tasks", "Patients", and "Logout". Below the navigation bar, there is a "Completed tasks" section listing various events and their dates. On the left side of the screen, there is a sidebar titled "Events" which lists several action items for different patients and dates. The main content area displays a "Paper Stack" for a specific event. The stack shows the date "4 Saturday August 2012" and the text "Patient 3 - 3 manually took 1 pill of the 1st dose." Below this, there is a list of "Steps to be completed:":

- Patient 3 - 3 was correctly identified by the system. (Icon: person)
- Pill 1 of 3 was taken correctly (Icon: red circle with 1)
- Pill 2 of 3 was taken correctly (Icon: red circle with 2)
- Pill 3 of 3 produced errors. Please review. (Icon: red circle with 3)
- Contact Patient 3 - 3 to discuss incident. (Icon: envelope)
- Close event.

At the bottom of the paper stack, there is a note: "Contact Patient 5 - 5 to discuss incident." There is also a "Settings" link at the bottom right of the main content area.

# Events Section - Using the Paper Stack

Selecting the Add Note icon will display a text box where notes may be added, such as comments received when contacting a patient about the action item.



Adding a note about an action item will move it to the PENDING category on the left.

De-identified video and images will be viewable for each medication administration event. Click the icon for each pill to expand and view.

# Events Section - Completed Tasks


Completed tasks list all tasks that have been completed by the user in the dashboard. These completed tasks may be sorted by:

- 1) Event Type, by selecting the EVENTS tab;
- 2) Task Type, by selecting the TASKS tab;
- 3) Patient, by selecting the PATIENTS tab; and
- 4) Date, by selecting the DATE tab (default).

The screenshot shows the AIUREX TECHNOLOGIES dashboard interface. At the top, there is a header bar with the text "Pharmaceutical 2012 Trial | Site NY-1 | Med A | Group 1" and "Pt: Dr. John Smith, PhD". To the right of the header are icons for a magnifying glass, a person, and a gear, followed by the text "AKUREX TECHNOLOGIES" and "Hello, Coordinator 1. Logout".

The main area of the dashboard has tabs for "Events", "Patients", and "Tasks". The "Events" tab is currently selected, showing a list of events for different patients and dates. One event is highlighted: "Patient 5 - 5 missed 1 pill of the 1st dose." on Sunday, September 30, 2012.

On the right side of the screen, there is a sidebar titled "Completed tasks" which lists various completed actions for different dates. The tasks include: "Completed steps for 8 - 8 from 2012-10-03.", "Reviewed pill 3 for 3 - 3 from 2012-09-28.", "Contacted 8 - 8 about 2012-10-03.", "Added note for 8 - 8 about 2012-10-03.", "Reviewed face identification for 7 - 7 from 2012-09-30.", "Reviewed pill 3 for 5 - 5 from 2012-09-30.", "Contacted 8 - 8 about 2012-10-02.", and "Contacted 8 - 8 about 2012-10-01.". The entire "Completed tasks" sidebar is highlighted with an orange border.



Ai-CURE TECHNOLOGIES

- 1 About AiDashboard
- 2 Using the Events Section
- 3 Using the Patients Section
- 4 Using the Patient Chart Section
- 5 Analytics and Settings

# Patients Section

Actuating the PATIENTS selector at the top left of the page will bring you to the PATIENTS HOME PAGE. (You can return to the EVENTS HOME PAGE at any time by actuating the EVENTS selector at the top left of the page.) The PATIENTS page is similarly divided into three main panels:

## PATIENT FILTERS

allow you to filter your patient list by patient type.

This panel contains a navigation bar with 'Events' and 'Patients' tabs, where 'Patients' is selected. Below this is a sidebar titled 'All patients' with three filter options: 'Patients with missed doses', 'Patients with late doses', and 'Patients who are late for current dose'. At the bottom is a search bar labeled 'search patients by name' with a magnifying glass icon and a '+' button.

## PATIENT SUMMARY DATA

Displays a list of all your patients, 7-day medication adherence rates, device battery levels, network signals, and allows you to message the patient..

This panel displays a summary for 10 patients. It includes a table with columns for Patient ID, Name, Age, Gender, Adherence Performance (last 7 days), Battery level, Signal strength, and Contact information. The adherence table shows a weekly grid with green squares for most days and red 'skip' squares for some. The battery and signal icons show varying levels of charge and strength. A search bar at the bottom is labeled 'search patients by name'.

Patient ID	Name	Age	Gender	Adherence Performance (last 7 days)							Battery	Signal	Contact
				Fr	Sa	Su	Mo	Tu	We	Today			
AI-001	1 - 1	25	male	green	green	green	green	max	86%	full	72%		
AI-002	2 - 2	25	male	green	green	green	green	-	100%	full	95%		
AI-003	3 - 3	25	male	green	green	green	green	green	100%	full	88%		
AI-004	4 - 4	25	female	green	green	green	green	-	100%	full	87%		
AI-005	5 - 5	30	male	green	green	skip	green	green	86%	full	92%		
AI-006	6 - 6	60	female	green	green	green	green	green	100%	full	33%		

## COMPLETED TASKS

log tasks that have been completed in the dashboard.

This panel lists completed tasks grouped by date. It includes a header with tabs for 'Events', 'Tasks', 'Patients', and 'Date'. The tasks are listed with checkmarks and descriptions. The dates shown are October 4, 2012, October 2, 2012, and October 1, 2012.

Completed tasks	
Thursday	Oct 4 2012
✓ Completed steps for 8 - 8 from 2012-10-03.	
✓ Reviewed pill 3 for 3 - 3 from 2012-09-28.	
✓ Contacted 8 - 8 about 2012-10-03.	
✓ Added note for 8 - 8 about 2012-10-03.	
✓ Reviewed face identification for 7 - 7 from 2012-09-30.	
✓ Reviewed pill 3 for 5 - 5 from 2012-09-30.	
Tuesday	Oct 2 2012
✓ Contacted 8 - 8 about 2012-10-02.	
Monday	Oct 1 2012
✓ Contacted 8 - 8 about 2012-10-01.	

# Patients Section – Patient Filters

Filter through the PATIENT SUMMARY DATA PANEL by using any of the preset filters to display:

- 1) ALL PATIENTS;
- 2) PATIENTS WITH MISSED DOSES;
- 3) PATIENTS WITH LATE DOSES; and
- 4) PATIENTS WHO ARE LATE FOR CURRENT DOSE.

Press the “+” button at the bottom left to customize your filter by GENDER, EVENT TYPE, and TO and FROM dates. You may also search for a particular patient by patient identifier.

The screenshot shows the Pharmaceutical 2012 Trial software interface. At the top, there is a navigation bar with the text "Pharmaceutical 2012 Trial | Site NY-1 | Med A | Group 1" and "Pt: Dr. John Smith, PhD". To the right of the navigation bar are icons for a magnifying glass, a double arrow, and a power button. The main header reads "10 patients." Below the header is a table with columns: Patient ID, Name, Age, Gender, Adherence Performance (last 7 days), Battery, Signal, and Contact. The table lists six patients (AI-001 to AI-006) with their respective details. On the far left, there is a sidebar with three buttons: "All patients" (highlighted with an orange border), "Patients with missed doses", and "Patients with late doses". At the bottom left of the main area is a small orange square containing a white plus sign (+). At the bottom center is a search bar with the placeholder "search patients by name". On the far right, there is a sidebar titled "Completed tasks" which lists various tasks with checkmarks and dates (e.g., Oct 4 2012, Oct 2 2012, Oct 1 2012).

Patient ID	Name	Age	Gender	Adherence Performance (last 7 days)	Battery	Signal	Contact
AI-001	1 - 1	25	male				
AI-002	2 - 2	25	male				
AI-003	3 - 3	25	male				
AI-004	4 - 4	25	female				
AI-005	5 - 5	30	male				
AI-006	6 - 6	60	female				

# Patients Section – Patient List


The PATIENT SUMMARY DATA displays:

- 1) Patient ID;
- 2) ADHERENCE PERFORMANCE for the last seven (7) days.
- 3) BATTERY LEVEL
- 4) DEVICE SIGNAL STRENGTH
- 5) Ability to MESSAGE the device

Select a PATIENT from the list to go directly to a detailed PATIENT REPORT.

The screenshot shows a software interface for managing patients in a clinical trial. At the top, there's a header with the trial information: "Pharmaceutical 2012 Trial | Site NY-1 | Med A | Group 1" and "PI: Dr. John Smith, PhD". To the right, it says "Hello, Coordinator 1." and "Logout". Below the header, there are tabs for "Events", "Patients", and "Completed tasks". The "Patients" tab is selected, showing a list of 10 patients. Each patient entry includes their ID, name, age, gender, adherence performance (a bar chart for the last 7 days), battery level, signal strength, and a message icon. Patient AI-002 is highlighted with an orange box. The "Completed tasks" section on the right lists various interactions with the patient, such as completed steps, reviews, and contacts, along with dates like October 4, 2012, and October 2, 2012.

Patient ID	Name	Age	Gender	Adherence Performance (last 7 days)	Battery	Signal	Contact
AI-001	1 - 1	25	male				
AI-002	2 - 2	25	male				
AI-003	3 - 3	25	male				
AI-004	4 - 4	25	female				
AI-005	5 - 5	30	male				
AI-006	6 - 6	60	female				



Ai-CURE TECHNOLOGIES

- 1 About AiDashboard
- 2 Using the Events Section
- 3 Using the Patients Section
- 4 Using the Patient Chart Section
- 5 Analytics and Settings

# Patient Chart

The PATIENT REPORT includes selection tabs about the patient SUMMARY, COMPLIANCE, NOTES, MESSAGES and SETTINGS. The following information about a particular patient is displayed:

## SUMMARY

includes PATIENT INFO and TRIAL RELATED data.

## COMPLIANCE

provides a DE-IDENTIFIED VIDEO AUDIT trail of each medication administration event.

## NOTE

displays any NOTES entered into the system about this patient.

## MESSAGE

displays any all communications to a patient's device.

## SETTINGS

allow changes to the patient late time.

The screenshot shows the AI-CURE TECHNOLOGIES Patient Chart interface. At the top, there is a navigation bar with tabs for "Events", "Patients", and "Date". On the left, a sidebar lists categories: "All patients", "Patients with missed doses", "Patients with late doses", and "Patients who are late for current dose". The main content area is titled "Patient 1 - 1". It features a header with five tabs: "Summary" (highlighted with an orange border), "Compliance", "Note", "Message", and "Settings". Below the tabs is a "Patient Information" section containing a placeholder image and the following details:

- Patient ID: AI-001
- Age, Gender: 36, male
- Today: 2012-08-10, C1D18 (18th day in 1st cycle)
- Treatment Start Date: 2012-07-24
- Treatment Stop Date: 2012-08-31

Under the "Current Trial" section, the following information is listed:

- Med A : 10mg/pill, 3 pills daily
- 1st dose : 10:00:00(reminder), 3 pills
- 1st cycle : 2012-07-24 ~ 2012-08-31
- Total # of "Miss Dose" : 1
- Total # of "Skip Dose" : 0
- Total # of "flagged" : 0
- Average minutes of "Late Dose" : 75

At the bottom of the main content area, there is a search bar with the placeholder "search patients by name" and a "+" button. To the right of the main content area, a vertical sidebar displays a log of events for the patient:

Day	Date	Event Details
Friday	Aug 10 2012	Completed steps for 8 - 8 from 2012-08-10. Added note for 8 - 8 about 2012-08-10. Completed steps for 8 - 8 from 2012-08-08. Completed steps for 7 - 7 from 2012-08-06. Contacted 7 - 7 about 2012-08-06. Reviewed pill3 for 7 - 7 from 2012-08-06. Reviewed face identification for 7 - 7 from 2012-08-06. Added note for 8 - 8 about 2012-08-10. Contacted 8 - 8 about 2012-08-10. Contacted 1 - 1 about 2012-08-09. Added note for 1 - 1 about 2012-08-09. Added note for 1 - 1 about 2012-08-09.
Wednesday	Aug 8 2012	Contacted 8 - 8 about 2012-08-08.
Tuesday	Aug 7 2012	Contacted 8 - 8 about 2012-08-07.

At the very bottom right, there is a "Settings" link.

# Patient Chart – Compliance

In Patient Compliance, you may visually review and flag every patient medication administration through a full video audit trail. Each of the icons below represents a video for each pill taken for each day the patient is on the trial. Below the video audit trail is a daily account of the patient's behavior.

-  denotes watched video
-  denotes flagged video
-  denotes video flagged by super user

Pharmaceutical 2012 Trial | Site NY-1 | Med A PI: Dr. John Smith, PhD

AiCURE TECHNOLOGIES

Hello, Coordinator 1. Logout

Events Patients

All patients

Patients with missed doses

Patients with late doses

Patients who are late for current dose

Patient 30105

Summary Compliance Note Message Settings

Print

Compliance Information - Video Audit

Oct	27	28	29	30	31	Nov	01	02	03	Today	05	06	07	08	09	10
					skip											

Medication: Med A


Date	Dose	Event	Flag
2012-11-04	1	Normal dose taken using AiView	
2012-11-03	1	Normal dose taken using AiView	
2012-11-02	1	Normal dose taken using AiView	

+ search patients by ID

Sunday Nov 4 2012

- ✓ Completed steps for patient 30105 from 2012-10-31.
- ✓ Contacted patient 30105 about 2012-10-31.
- ✓ Contacted patient 30103 about 2012-10-29.
- ✓ Added note for patient 30105 about 2012-10-31.
- ✓ Reviewed pill3 for patient 30105 from 2012-10-31.

Settings



Ai-CURE TECHNOLOGIES

- 1 About AiDashboard
- 2 Using the Events Section
- 3 Using the Patients Section
- 4 Using the Patient Chart Section
- 5 Analytics and Settings

# Analytics

Through the dashboard's ANALYTICS section, you may access aggregated medication adherence data for the entire patient population within a site. This ANALYTICS section may be accessed by clicking on the orange arrows at the top header on each page to gain insight into patient performance and trends.

The screenshot displays two views of the AI-CURE TECHNOLOGIES dashboard. The top view shows the main interface with a navigation bar, event and patient tabs, and a central message about a patient missing doses. The bottom view shows the 'Overall Patient Adherence' section with a pie chart illustrating usage statistics. An orange double-headed arrow points between the two sections, indicating they are interconnected.

**Top Dashboard View:**

- Header: Pharmaceutical 2012 Trial | Site NY | Med A | Group 1 PI: John Smith, PhD
- Header Buttons: Up and Down arrows (highlighted with an orange box), AI-CURE TECHNOLOGIES logo, Hello, Coordinator 1, Logout
- Tab Bar: Events, Patients, Tasks (highlighted), Date
- Completed tasks: Friday, Aug 10 2012
  - ✓ Completed steps for 8 - 8 from 2012-08-10.
  - ✓ Added note for 8 - 8 about 2012-08-10.
  - ✓ Completed steps for 8 - 8 from 2012-08-08.

**Bottom Dashboard View:**

- Header: Pharmaceutical 2012 Trial | Site NY Med A Group 1 PI: John Smith, PhD
- Header Buttons: Up and Down arrows (highlighted with an orange box), AI-CURE TECHNOLOGIES logo, Hello, Coordinator 1, Logout
- Section: Overall Patient Adherence
- Section: Overall Patient Usage (Pie Chart)
  - Successful, 97.18%
  - Missed, 2.26%
  - Skipped, 0.56%
- Legend: Successful (green), Skipped (orange), Missed (red)

# Settings

Within the USER SETTINGS, users may adjust settings related to patient alerts. Contact methods and alert thresholds may also be changed to ensure notifications are sent according to individual user preferences. Accountability events – if a patient misses a certain number of doses within a specific period of time - may also be enabled. These in turn will generate ACTION ITEMS

The ADD PATIENTS tab may be used to enter patients, who are under your supervision, into the trial.

The screenshot shows the 'User Settings' page of the AI-CURE TECHNOLOGIES application. At the top, there is a navigation bar with the text 'Pharmaceutical 2012 Trial | Site NY | Med A | Group 1' and a dropdown menu, followed by 'PI: John Smith, PhD'. To the right of the PI information is the AI-CURE TECHNOLOGIES logo and a 'Hello, Coordinator 1.' message with a 'Logout' link. Below the navigation bar, there are two tabs: 'Add Patients' (which is currently selected) and 'User Settings'. The main content area contains the following configuration options:

- Contact Information:** Fields for Phone (555-555-5555), Email (email@email.com), and Title (Mr.).
- Preferred Contact Method:** Radio buttons for None, Text (selected), and Email.
- Preferred Contact Time:** A dropdown menu set to 12:00.
- Alert Thresholds:** Checkboxes for 'Contact me if patient adherence drops below [60%]' over the 'whole trial' and 'Contact me if patient misses [3]' days in a row.
- Event Generation:** Checkboxes for 'Generate Events For': 'If patient misses [3]' doses in a row and 'If patient misses [6]' doses in the past '[14]' days.
- Grouping:** A checkbox for 'Group similar event types in the Events module'.
- First Loading Page:** Radio buttons for 'Event Page' (selected) and 'Patient List' (with a note '(Event page is the default loading page)').

At the bottom of the form is a 'Save' button. Below the form, a link '« Back to dashboard' is visible.



# AiDashboard™ Summary

Designed to enable optimal medication adherence

Access real-time medication adherence data

De-identified video audit of each medication administration

Intervene with patients before they deviate from protocol

URL: [www.myclinicaltrialdashboard.com](http://www.myclinicaltrialdashboard.com)

Call Technical Support: +1 (800) 675-4960 [7am to 7pm PST]

Email Technical Support: [support@aicuretechnologies.com](mailto:support@aicuretechnologies.com)

A full user manual is available through the help button at the top right of each page.

