




AiDashboard™

Superuser Guide

AI-CURE TECHNOLOGIES



AI-CURE TECHNOLOGIES

URL: www.myclinicaltrialdashboard.com


Call Technical Support: +1 (800) 675-4960 [7am to 7pm PST]

Email Technical Support: support@aicuretechnologies.com

Name:

Username:

Password:



Ai-CURE TECHNOLOGIES

1 About AiDashboard

2 Using the Patients Section


3 Using the Patient Chart Section

4 Analytics


About the AiDashboard™

The AiDashboard™ allows you to track study participant performance. The dashboard provides real time adherence data for study participants, indicates which participants may need help, and suggests where study coordinators can best spend their time. Complete audit trails for each participant are also available.


Quick Reference



Participant & Study Coordinator Daily Use



www.myclinicaltrialdashboard.com



Ai-CURE TECHNOLOGIES

1 About AiDashboard

2 Using the Patients Section

3 Using the Patient Chart Section


4 Analytics

AiDashboard Login

Each super user will be provided with a unique username and password to access the secure AiDashboard. To logon, please go to:

www.myclinicaltrialdashboard.com

If you forget your password, please click on the orange “Forgot your password?” link, enter your email address and username on file, and a new password will be sent to your address.



Patients Section

When you first login, you will arrive at the PATIENTS HOME PAGE, which shows adherence data for all patients in the trial. This PATIENTS HOME PAGE shows summary data on the left, with the full patient list broken down by site in the middle section. In these you will find the following information:

SUMMARY DATA

Gives a quick look at overall trial data and rankings.

Pharmaceutical 2012 Trial | Med A | All Patients ▾

Overall Patient Statistics All patients in trial	
Overall Adherence:	89.59%
Total # of patients:	100
Number of days into trial:	10
Highest performing site:	10-IL-2 (96.91%)
Lowest performing site:	06-CA-3 (61.86%)

 Print this page

Print







PATIENT SUMMARY DATA

Displays a list of all patients in the trial, broken down by site (clicking on the “+” on the right-hand side of each site expands its data). 7-day medication adherence rates and device battery levels are also shown.



AI-CURE TECHNOLOGIES

Hello, Super User. Logout

01-NY-1		PI: Dr. John Smith, PhD		Overall Compliance: 94.84%	
Coordinator 1 - Overall Compliance: 94.85%		Videos reviewed: 0		Coordinator 2 - Overall Compliance: 100%	
Videos reviewed: 0		Videos reviewed: 0		Videos reviewed: 0	
Patient ID	# Days into trial	Overall Compliance	Managed By	Adherence Performance (last 7 days)	Battery
30101	10	90%	Coordinator 1	 miss	 59%
30102	10	100%	Coordinator 1	 late	 85%
30103	10	100%	Coordinator 1, Coordinator 2	 m	 24%
30104	10	100%	Coordinator 1	 -	 94%
30105	10	90%	Coordinator 1	 skip	 26%
30106	10	100%	Coordinator 1		 56%
30107	10	100%	Coordinator 1	 m	 86%
30108	10	67%	Coordinator 1	 miss miss miss late	 68%
30109	10	100%	Coordinator 1		 75%
30110	10	100%	Coordinator 1, Coordinator 2		 77%

02-NY-2

PI: Dr. John Smith, PhD

Overall Compliance: 95.87% +

Patients Section – Summary Data

Left column “score card” shows summary data about the trial:

- 1) Overall Adherence;
- 2) # of Patients;
- 3) # of Days Into the Trial;
- 4) Highest Performing Site;
- 5) Lowest Performing Site

Pharmaceutical 2012 Trial | Med A | All Patients

AI-CURE TECHNOLOGIES Hello, Super User. Logout

Overall Patient Statistics
All patients in trial

Overall Adherence: 89.59%

Total # of patients: 100





Number of days into trial: 10

Highest performing site: 10-IL-2 (96.91%)

Lowest performing site: 06-CA-3 (61.86%)

Print this page

Print

01-NY-1		PI: Dr. John Smith, PhD		Overall Compliance: 94.84%	
				Videos reviewed: 0	Videos reviewed: 0
Patient ID	# Days into trial	Overall Compliance	Managed By	Adherence Performance (last 7 days)	
30101	10	90%	Coordinator 1		59%
30102	10	100%	Coordinator 1		85%
30103	10	100%	Coordinator 1, Coordinator 2		24%
30104	10	100%	Coordinator 1		34%
30105	10	90%	Coordinator 1		26%
30106	10	100%	Coordinator 1		56%
30107	10	100%	Coordinator 1		85%
30108	10	67%	Coordinator 1		68%
30109	10	100%	Coordinator 1		75%
30110	10	100%	Coordinator 1, Coordinator 2		77%

02-NY-2 PI: Dr. John Smith, PhD Overall Compliance: 95.87%

Patients Section – Patient List

The PATIENT SUMMARY DATA displays:

- 1) Patient ID;
- 2) Number of DAYS into the trial;
- 3) OVERALL COMPLIANCE;
- 4) Who the patient is MANAGED BY;
- 5) ADHERENCE PERFORMANCE for the last seven (7) days;
- 6) BATTERY LEVEL

Pharmaceutical 2012 Trial | Med A | All Patients

AI-CURE TECHNOLOGIES

Hello, Super User. Logout

Overall Patient Statistics
All patients in trial

Overall Adherence: 89.59%
Total # of patients: 100
Number of days into trial: 10
Highest performing site: 10-IL-2 (96.91%)
Lowest performing site: 06-CA-3 (61.86%)

Print this page

Patient ID	# Days into trial	Overall Compliance	Managed By	Adherence Performance (last 7 days)	Battery
30101	10	90%	Coordinator 1		
30102	10	100%	Coordinator 1		
30103	10	100%	Coordinator 1, Coordinator 2		
30104	10	100%	Coordinator 1		
30105	10	90%	Coordinator 1		
30106	10	100%	Coordinator 1		
30107	10	100%	Coordinator 1		
30108	10	67%	Coordinator 1		
30109	10	100%	Coordinator 1		
30110	10	100%	Coordinator 1, Coordinator 2		


Coordinator 1 - Overall Compliance: 94.85%
Coordinator 2 - Overall Compliance: 100%
Videos reviewed: 0
Videos reviewed: 0

Overall Compliance: 94.84%

02-NY-2

PI: Dr. John Smith, PhD

Overall Compliance: 95.87%



Ai-CURE TECHNOLOGIES

1 About AiDashboard

2 Using the Events Section

3 Using the Patients Section

4 Using the Patient Chart Section

5 Analytics

Patient Chart

The PATIENT REPORT includes selection tabs about the patient SUMMARY, COMPLIANCE, and NOTES. The following information about a particular patient is displayed:

SUMMARY

includes PATIENT INFO and TRIAL RELATED data.

COMPLIANCE

provides a DE-IDENTIFIED VIDEO AUDIT trail of each medication administration event.

NOTE

displays any NOTES about flagged videos entered into the system about this patient.

Pharmaceutical 2012 Trial | Site NY-1 | Med A PI: Dr. John Smith, PhD

AiCURE TECHNOLOGIES Hello, Super User. Logout

Overall Patient Statistics
All patients in trial

Overall Adherence:	89.59%
Total # of patients:	100
Number of days into trial:	10
Highest performing site:	Pharmaceutical 2012 Trial Site IL-2 (96.91%)
Lowest performing site:	Pharmaceutical 2012 Trial Site CA-3 (61.86%)

Patient 30104

Summary **Compliance** **Note**

Patient Information

	Patient ID: 30104
	Today: 2012-11-16, C1D10 (10th day in 1st cycle)
	Treatment Start Date: 2012-11-07
	Treatment Stop Date: 2012-12-31

Current Trial

Med A	: 10mg/pill, 3 pills daily
1st dose	: 10:00:00(reminder), 3 pills
1st cycle	: 2012-11-07 ~ 2012-12-31
Total # of "Miss Dose"	: 0
Total # of "Skip Dose"	: 0
Total # of "flagged"	: 2
Average minutes of "Late Dose"	: 117

« Back to patient list

Patient Chart – Compliance

In Patient Compliance, you may visually review and flag every patient medication administration through a full video audit trail. You can also see which videos have been reviewed by a study coordinator. Each of the icons below represents a video for each pill taken for each day the patient is on the trial. Below the video audit trail is a daily account of the patient's behavior.


The screenshot displays the 'Patient 30104' chart within the 'Pharmaceutical 2012 Trial | Site NY-1 | Med A' section. The top navigation bar includes 'PI: Dr. John Smith, PhD', 'A-CURE TECHNOLOGIES', and a 'Logout' link. The left sidebar shows 'Overall Patient Statistics' for all patients in trial, including Overall Adherence (89.59%), Total # of patients (100), Number of days into trial (10), Highest performing site (Pharmaceutical 2012 Trial | Site IL-2 (96.91%)), and Lowest performing site (Pharmaceutical 2012 Trial | Site CA-3 (61.86%)).

The main chart area for Patient 30104 features a title 'Patient 30104' and tabs for 'Summary', 'Compliance', and 'Note'. The 'Compliance' tab is active, showing a 'Compliance Information - Video Audit' grid for November. The grid has columns for dates from Nov 08 to Nov 30. Rows represent individual patient pills, with icons indicating whether they were watched or flagged by a coordinator or super user. A legend at the top right defines the icons: blue square for 'Watched by coordinator', green square for 'Watched by super user', yellow square for 'Watched by both', red triangle for 'Flagged by coordinator', and orange triangle for 'Flagged by super user'.

The 'Medication: Med A' table below the grid lists medication events:

Date	Dose	Event	Flag
2012-11-16	1	Awaiting data from AiView	
2012-11-15	1	194 minutes late but dose taken	
2012-11-14	1	Normal dose taken using AiView	
2012-11-13	1	Normal dose taken using AiView	
2012-11-12	1	Normal dose taken using AiView	
2012-11-11	1	Normal dose taken using AiView	🚩

At the bottom of the chart area is a link '« Back to patient list'.



Ai-CURE TECHNOLOGIES

1 About AiDashboard

2 Using the Events Section

3 Using the Patients Section

4 Using the Patient Chart Section

5 Analytics

Analytics

Through the AiDashboard's ANALYTICS section, you may access aggregated medication adherence data for the entire patient population, or compare adherence rates by site or by coordinator. This ANALYTICS section may be accessed by clicking on the orange arrows at the top header on each page to gain insight into patient performance and trends.

The screenshot illustrates the transition from the main dashboard to the Analytics section. At the top, an orange double-headed arrow highlights the header area. Below it, the main dashboard displays overall patient statistics (Overall Adherence: 89.59%) and site-specific details for two sites: 01-NY-1 and 02-NY-2. Site 01-NY-1 has an Overall Compliance of 94.84%, while Site 02-NY-2 has 95.87%. The Analytics section below provides a detailed breakdown of patient adherence, including a pie chart showing usage across three categories: Successful (89.59%), Skipped (6.29%), and Missed (4.12%). The overall statistics summary on the right reiterates the 89.59% adherence rate and includes the total number of patients (100) and the number of days into trial (10).

Pharmaceutical 2012 Trial | Med A | All Patients

AI•CURE•TECHNOLOGIES

Hello, Super User. Logout

Overall Patient Statistics
All patients in trial
Overall Adherence: 89.59%

01-NY-1 PI: Dr. John Smith, PhD Overall Compliance: 94.84% +

02-NY-2 PI: Dr. John Smith, PhD Overall Compliance: 95.87% +

Overall Patient Adherence ✓

Patient Analysis by Site

Patient Analysis by Coordinator

Overall Patient Usage Across All Sites

Successful, 89.59% Missed, 4.12% Skipped, 6.29%

Overall Statistics

Overall Adherence:	89.59%
Total # of patients:	100
Number of days into trial:	10

Analytics – By Site

Selecting Patient Analysis by Site in the left-hand column will show comparative data across all sites. The right hand-column lists the sites in descending order of compliance rate.



Analytics – By Coordinator

Selecting Patient Analysis by Coordinator in the left-hand column will show comparative data of all study coordinators across all sites. The right hand-column lists the study coordinators in descending order of compliance rate of the patients under their care.

Pharmaceutical 2012 Trial | Med A | All Patients

AI-CURE TECHNOLOGIES

Hello, Super User. Logout

Overall Patient Adherence

Patient Analysis by Site

Patient Analysis by Coordinator ✓

01-NY-1
Comparison of Coordinators

Type of Doses	Coordinator 1 (%)	Coordinator 2 (%)
Successful	94.84%	100%
Skipped	4.16%	0%
Missed	0%	0%

02-NY-2
Comparison of Coordinators

Type of Doses	Coordinator 3 (%)	Coordinator 4 (%)
Successful	96.55%	94.87%
Skipped	2.72%	3.13%
Missed	0%	0%

03-NY-3
Comparison of Coordinators

Type of Doses	Coordinator 9 (%)	Coordinator 10 (%)	Coordinator 11 (%)
Successful	93.75%	86.2%	85.71%
Skipped	5.17%	8.62%	8.16%
Missed	1.08%	5.11%	6.12%

04-CA-1
Coordinator 5

Type of Doses	Percentage
Successful	91.75%
Skipped	4.12%
Missed	4.12%

05-CA-2
Comparison of Coordinators

Type of Doses	Coordinator 6 (%)	Coordinator 7 (%)
Successful	97.8%	95.88%
Skipped	1.08%	3.11%
Missed	0%	0%

06-CA-3
Coordinator 8

Type of Doses	Coordinator 8 (%)
Successful	97.8%
Skipped	1.08%
Missed	0%

Coordinators listed in order of compliance

Coordinator (Site)	Compliance Rate
Coordinator 2 (01-NY-1)	100%
Coordinator 6 (05-CA-2)	98.28%
Coordinator 18 (10-IL-2)	97.96%
Coordinator 3 (02-NY-2)	96.55%
Coordinator 17 (10-IL-2)	95.83%
Coordinator 13 (07-CA-4)	95.4%
Coordinator 12 (07-CA-4)	95%
Coordinator 4 (02-NY-2)	94.87%
Coordinator 1 (01-NY-1)	94.85%
Coordinator 15 (08-CA-5)	93.88%
Coordinator 9 (03-NY-3)	93.75%
Coordinator 7 (05-CA-2)	93.22%
Coordinator 16 (09-IL-1)	92.78%
Coordinator 5 (04-CA-1)	91.75%
Coordinator 10	

AiDashboard™ Summary

Designed to enable optimal medication adherence

Access real-time medication adherence data

De-identified video audit of each medication administration

Study coordinators intervene with patients before they deviate from protocol

URL: www.myclinicaltrialdashboard.com

(supported browsers):

Internet Explorer	Firefox*	Chrome*	Safari*	Opera
8+	4+	any	5+	12+

* Recommended

Call Technical Support: +1 (800) 675-4960 [7am to 7pm PST]

Email Technical Support: support@aicuretechnologies.com

A full user manual is available through the help button at the top right of each page.

